

Magen R4130与其它品牌试剂的对比

Company	样品用量	核酸浓度	Unit	260/280	260/230	Yield	µg
Magen R4130	肝(50mg)	0.362	µg/µl	2.14	2.12	173.8	
	肝(50mg)	0.335	µg/µl	2.13	2.16	160.7	
Qiagen	肝(50mg)	0.333	µg/µl	2.13	1.66	159.9	
	肝(50mg)	0.345	µg/µl	2.12	1.65	165.7	
Invitrogen	肝(50mg)	0.337	µg/µl	2.13	1.96	162.0	
	肝(50mg)	0.360	µg/µl	2.13	1.75	172.8	
Omega	肝(50mg)	0.328	µg/µl	2.14	1.73	157.3	
	肝(50mg)	0.360	µg/µl	2.15	1.71	172.8	
Company	样品用量	核酸浓度	Unit	260/280	260/230	Yield	µg
Magen R4130	心(50mg)	0.157	µg/µl	2.12	1.34	15.7	
	心(50mg)	0.141	µg/µl	2.13	1.52	14.1	
Qiagen	心(50mg)	0.148	µg/µl	2.13	1.17	14.8	
	心(50mg)	0.149	µg/µl	2.11	0.95	14.9	
Invitrogen	心(50mg)	0.133	µg/µl	2.12	1.43	13.3	
	心(50mg)	0.153	µg/µl	2.11	1.39	15.3	
Omega	心(50mg)	0.154	µg/µl	2.06	0.73	15.4	
	心(50mg)	0.131	µg/µl	2.1	1.15	13.1	